


RENA - MULTICART[®]

Reusable Bicarbonate cartridge for Haemodialysis

Available for first time



*

For use with Fresenius NG Machines



*

For use with Nipro Surdial Machines

Just open the removable cap, rinse, fill the Bicarbonate Granules & reuse the cartridge.
Advantages of Bicarbonate Cartridge:

1. Delivers Bicarbonate dosage exactly as prescribed unlike liquid Bicarbonate which deteriorates with time and converts partially to carbonate thereby delivering a much lower dosage of Bicarbonate than prescribed.
2. Eliminates risk of Bacterial contamination as Bicarbonate Solution is prepared online.
3. Hygienic, easy & convenient method of Bicarbonate Conc. Supply.

Available in two sizes - 550/650 gms granulated, free flow sodium bicarbonate powder.

Manufactured and Marketed by:
ATLANTIC PHARMACEUTICALS

101, Business Point 349, Western Express Highway, Andheri East, Mumbai - 400 069
Ph +91 22 61594343/26820950 : Email - atlanticpharma@vsnl.com